# SuperSeton

Smooth, simple, secure.

CE 1639

## USER MANUAL

EN	User manual	3	FI	Käyttöopas	36
DE	Bedienungsanleitung	6	PL	Instrukcja użytkownika	39
FR	Manuel utilisateur	9	HU	Használati útmutató	42
NL	Gebruikshandleiding	12	CS	Uživatelská příručka	45
ES	Manual de usuario	15	BG	Ръководство за потребителя	48
IT	Manuale d'uso	18	HR	Korisnički priručnik	51
PT	Manual do utilizador	21	RO	Manual de utilizare	54
GR	Εγχειρίδιο χρήστη	24	LV	Lietotāja rokasgrāmata	57
DA	Brugervejledning	27	LT	Naudotojo vadovas	60
SV	Användarhandbok	30	SL	Uporabniški priročnik	63
NO	Håndbok	33	SK	Používateľská príručka	66

Figure 1


Figure 2


Figure 3

Figure 4



Figure 5

# USER MANUAL

## 1 INDICATION FOR USE

The intended use of the SuperSeton is the drainage of perianal fistulae.

## 2 CONTRAINDICATIONS

- Do not use the SuperSeton in patients with an allergy to materials contained in the SuperSeton (medical grade silicon, medical grade polypropylene).

## 3 WARNINGS

- This product is not intended for other procedures than the drainage of anal fistulas.
- Do not use if the product has passed the expiration date printed on the product label.
- Inspect the sterile packaging to ensure it is intact. Do not use if the sterile packaging has been opened or damaged.
- This product is for single-use only, do not re-sterilize or re-use any part of this product.
- Do not attempt to alter this product in any way. Doing so may endanger the patient and / or user.
- Shorten the tube only at the end to which no insert is attached.
- Do not remove the insert from the tube.
- In order to close the SuperSeton, grasp the tube end with the applier and push the insert into the tube end. Do not pull on the tube during assembly.
- Align the insert and Applier before applying force to prevent bending and or buckling the insert.
- For secure closure of the SuperSeton, push the tube over the insert until both ends of the tube are pushing directly against each other and the insert is fully covered by the tube.
- SuperSeton may open under excessive tensile strain. This may result in a

prematurely opened SuperSeton which has to be replaced.

- Advise the patient to consult a medical practitioner if he or she has any unusual symptoms. Also when premature opening of the SuperSeton occurs.
- Give notice to the patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the patient is established.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the physicians or medical professional is established.
- The SuperSeton is intended for short-term use for a period no longer than 29 days.

## 4 PRECAUTIONS

- Read all instructions prior to use.
- Only physicians or medical professionals qualified in the appropriate surgical techniques should use this device.
- This product represents a potential biohazard after use. Handle and dispose of the product as required by hospital policy and applicable laws.

## 5 INCLUDED IN PACKAGE

See **Figure 1**.

- Applier
- Insert
- Drainage tube

## 6 INSTRUCTIONS FOR USE

- Check the package for damage, do not use the SuperSeton if the package is damaged.
- Make sure there is a spare SuperSeton pouch available.
- Open the package in a sterile manner and take the product out of the package.
- Introduce the tube end without the insert into the fistula tract with forceps.

5. Feed the tube end through the fistula, until the tip can be felt at the other end of the fistula tract.
6. Grasp the end of the tube and pull out of the fistula, so that both ends of the tube can still be felt; one at the entrance and one at the exit of the fistula.
7. Place the end of the tube without the insert into the applier, with the arrow on the applier pointing in the direction of the fistula tract (see **Figure 2**).
8. Shorten the tube to fit the patient (see **Figure 3**).
  - **Warning:** Shorten the tube at the end without the insert. Make sure the length of the tube is such that it can be closed without resulting in a permanent tension on the resulting connection.
  - When the tube is left too long the resulting smooth loop has a higher risk of being pulled in daily use. The optimal length results in the seton being close to the body without tension on the connection.
9. Push the insert into the tube to close the seton into a smooth loop. Push the insert in the direction of the arrow into the tube enclosed in the Applier (see **Figure 4**).
  - Do not pull on the tube or the insert.
  - For a secure closure of the drainage catheter the open end of the tube must fully cover the insert (see **Figure 5**).
  - The SuperSeton is intended for short-term use for a period no longer than 29 days.
10. When removing the SuperSeton do not cut on or close to the insert area.

Please review our general terms and conditions on [www.superseton.com](http://www.superseton.com)

## 7 SYMBOLS

Symbol	Description	Symbol	Description
	Manufacturer		Date of manufacture
	CE-mark and identification number of Notified Body.		Caution
	Use-by date		Consult instructions for use
	Batch number		Keep away from sunlight
	Catalogue number		Do not use if package is damaged
	Sterilized using irradiation		Do not re-use
	Do not re-sterilize		Temperature limit
	Single <i>sterile</i> barrier system		Medical device
	Double <i>sterile</i> barrier system		Unique device identifier
	Keep dry		

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**Super Seton B.V.**