

PRODUCT SPECIFICATION & DATA

PRODUCT:	Drain-Guard
REF:	8150

PRODUCT DIMENSIONS:	
Total Length of packed product:	230 mm
Product size	85 mm x 84 mm

PRODUCT COMPONENTS		
Component No.	Description	Material/Data
1	Backing paper	Paper + PE coated with Lina silicone system 8500
2	Primary carrier	Hydrocolloid and polyurethane film
3	Adhesive for cushion	Acrylic adhesive with thin paper carrier
4	Support cushion	Polyethylene foam
5	Adhesive for catheter friction	Acrylic adhesive with thin paper carrier
6	Siliconized paper protection	Paper + PE coated with Lina silicone system 8500
7	Non-woven / hydrocolloid lid	Polypropylene with acrylic adhesive combined with hydrocolloid

Packaging	1 pcs in sterile peel-pouch 150 x 230 mm 25 pcs in box 240 x 92 x 155 mm 4 x 25 pcs in export carton 384 x 252 x 160 mm
Sterilisation	Gamma sterilisation



