



## PRODUCT SPECIFICATION & DATA

PRODUCT:	Epi-Guard
REF:	8170

PRODUCT DIMENSIONS:	
Total Length of packed product:	165 mm
Product size	54 mm x 53 mm

PRODUCT COMPONENTS		
Component No.	Description	Material/Data
1	Backing paper	Paper + PE coated with Lina silicone system 8500
2	Primary carrier	Hydrocolloid and polyurethane film
3	Support cushion	Polyethylene foam with adhesive
4	Adhesive for catheter friction	Acrylic adhesive with PVC carrier
5	Siliconized paper protection	Paper + PE coated with Lina silicone system 8500
6	Non-woven / hydrocolloid lid	Polypropylene with acrylic adhesive combined with hydrocolloid

Packaging	1 pcs in sterile peel-pouch 100 x 165 mm 25 pcs in box 165 x 100 x 110 mm 4 x 25 pcs in export carton 342 x 210 x 118 mm
Sterilisation	Gamma sterilisation

