



LiNA Laparo Trainer

*Portable trainer for realistic training
in laparoscopic surgery*

- Quick and easy to transport and unfold
- Training at any time and anywhere
- Training with the same surgical instruments as used in the OR
- Eliminates the need for borrowing expensive equipment from the OR
- Lightweight aluminium suitcase
- Low investment – very cost-effective training

LiNA Laparo Trainer

The LiNA Laparo Trainer is developed to meet the growing demand for realistic laparoscopic training for young surgeons/students at an affordable price.

The LiNA Laparo Trainer is built in a portable lightweight aluminium suitcase - easy to move around or to take home and perform exercises or prepare for surgery next day. All you need apart from the LiNA Laparo Trainer is a table, a 220V socket and a few standard laparoscopic instruments. The LiNA Laparo Trainer is easily opened/unfolded and the flat screen is elevated to eye level in order to make your working position identical with the position you take in the operating room.

A digital video camera placed in the LiNA Laparo Trainer makes a perfect sharp image at a reasonable price. The training with the LiNA Laparo Trainer is very realistic since you are using the same surgical instruments as used in the operating room - a good way to get familiar with the laparoscopic instruments.



Training Kit

The LiNA Laparo Trainer comes with a training kit consisting of 8 practical training exercises. It is easy to make more exercise kits for further training sessions as they are made of materials available at all hospitals. It is, however, also possible to purchase additional training exercise kits. You will find up-dated information about the exercises at this webpage: www.laparotrainer.com

Ordering Information:

LiNA Laparo Trainer

REF.	Product	Packing
LT-1000	LiNA Laparo Trainer - Weight: 15,8kg. Dimensions: width 42cm x depth 24cm x height 67 cm.	1 pc.
LT-KITS	LiNA Laparo Trainer Kit - 8 training excercises.	1 pc.

LiNA distributor in your country:

