




**Sugita AVM Microclips** are made of cobalt-chromium-molybdenum-based alloy having superior properties, specified under ASTM F1058 G2(ISO5832-7).

Being made of a non-ferromagnetic alloy, they are safe for MRI currently in clinical use.

To prevent damage and to facilitate handling, the clip is individually seated in a heat-resistant silicone sponge and packaged in a metal container, which is fully autoclavable.


Clips are made to have a holding force within the predetermined range of 50-70g.



# Sugita A V M Micro Clip and Applier




*Clip Packaging*




# Sugita A V M Micro Clip

**MIZUHO®**


Catalog No.	Blade Length
07-932-02	2 mm
07-932-03	3 mm
07-932-04	4 mm
07-932-05	5 mm



*Clip can be easily picked up  
by an applier while holding  
the container in one hand.*



## Sugita A V M Micro Clip Applier



Catalog No.	Overall Length
07-935-10	21.2 cm

### Mizuho Medical Co., Ltd.

ICN building 5F 27-17,2-chome, Hongo, Bunkyo-ku,  
Tokyo, 113-0033 Japan  
Fax No. +81-3-3818-1705 Phone No. +81-3-3815-7101  
E-Mail : marketing@mizuhomedical.co.jp

### Mizuho America Inc.

133 Brimbal Avenue Beverly, MA 01915  
Phone No. 800-699-CLIP(2547) or 978-921-1718  
Fax No. 978-921-4003

Manufacturer

**MIZUHO IKAKOGYO CO., LTD.**  
Tokyo, Japan