



## PRODUCT SPECIFICATION & DATA

<b>PRODUCT:</b>	Pleura-Guard
<b>REF:</b>	8161

<b>PRODUCT DIMENSIONS:</b>	
<b>Total Length of packed product:</b>	230 mm
<b>Product size</b>	70 mm x 100 mm

<b>PRODUCT COMPONENTS</b>		
<b>Component No.</b>	<b>Description</b>	<b>Material/Data</b>
1	Backing paper	Paper + PE coated with Lina silicone system 8500
2	Primary carrier	Hydrocolloid and polyurethane film
3	Adhesive for cushion	Acrylic adhesive with PVC carrier
4	Support cushion	Polyethylene foam
5	Non-woven / hydrocolloid lid	Polypropylene with acrylic adhesive combined with hydrocolloid

<b>Packaging</b>	1 pcs in sterile peel-pouch 150 x 230 mm 25 pcs in box 298 x 189 x 189 mm 4 x 25 pcs in export carton 597 x 397 x 311 mm
<b>Sterilisation</b>	Gamma sterilisation

A blue ink handwritten signature, appearing to be "PV", is written over the "Approved by : PV" text.

