



SPIGGLE & THEIS
Medizintechnik

MENTIO

Medical ENT Instruments Overath

ENT-Single Use Instruments



25 Years

SPIGGLE & THEIS
1994 - 2019



MENTIO – ENT-Single Use Instruments

MENTIO


Disposable instruments - an economical and safe alternative

Single use instruments are an excellent solution for meeting both, the high demand for hygiene and safety for medical instruments and the economic requirements for efficiency and time.


The advantages of disposable instruments are self evident. In addition to the costs of validation, documentation, processing, and repair, disposable instruments eliminate liability risks. At the same time, the safety of you and your patients increases as there is no residual contamination.

Our MENTIO line is based on many years of experience in manufacturing and selling high-end ENT-instruments. Convince yourself and experience the high quality and precision of a new generation of disposable instruments.

The ideal solution for your practice and ward - **try ours!**



Single Use Instruments



- Disposable instruments for the highest quality demands
- Matt surface for reduced light reflection



- Save yourself reprocessing time
- Ready to use
- No preparation, repairs or logistics required



- No residual contamination

Nasal speculum according to Hartmann

Large

- blade length 3.2 cm
- total length 15.8 cm
- polyamide
- 25 pieces per box

5002002

Small

- blade length 2.6 cm
- total length 14.5 cm
- polyamide
- 25 pieces per box

5002001

matt surface for reduced light reflection

ergonomic grip

ribbed for better control




Single Use Instruments

Cerumen hook

- hook size 0.3 cm
- total length 13.5 cm
- carbon fibre filled nylon
- 100 pieces per box

5001007



Ear curette / Cotton carrier

Combination instrument with round loop and cotton wool carrier

- length cotton wool carrier 1.1 cm
- total length 15 cm
- carbon fibre filled nylon
- 100 pieces per box

5001008





Not suitable for use with otoscopes

Ear speculum according to Hartmann

- matt inner surface
- holding ring
- polypropylene
- 200 pieces per box

Art.No.	Type/Size	ID (mm)	OD (mm)
5001001	0	2.1	3.4
5001002	1	3.2	4.5
5001003	2	4.0	5.3
5001004	3	5.2	6.5
5001005	4	6.5	7.8



Acrylic Dispenser

- dispenser for disposable ear specula
- 5 compartments, 0.3 cm acrylic glass
- labeled inserts (size 0-4)

5001000



Spatula, birch wood

- 15 cm
- 100 pieces per box

98-11-402



© 2019 by SPIGGLE & THEIS Medizintechnik GmbH
Modifications and errors - including pictures -
reserved. Please note also that not all goods might
be for sale in your country.

FL_Mentio_E_01



Burghof 14 · 51491 Overath · Germany
Phone: +49 2206 90 81 - 0
Fax: +49 2206 90 81 - 13
info@spiggle-theis.com
www.spiggle-theis.com