

CINCH™

Tube Securement Device

9F - 30F


Comfort, First & Foremost

CINCH™ was designed with the most delicate skin in mind. Our hypoallergenic silicone adhesive provides a gentle yet secure attachment. The low-profile and flexible design of the CINCH™ makes it comfortable and non-disruptive, even when it's not in use.



Safe, Secure, Securement


CINCH™ promotes the health, safety, and comfort of patients by securely anchoring medical tubes while relieving stress and strain at insertion points.




CLEAN




APPLY




PLACE



SECURE



9F



12F



30F


CAUTION: If redness or irritation occurs, stop use and consult a healthcare professional. Do not over tighten. This document is not a substitute for the DFU (Directions for Use).
IMAGES ARE NOT TO SCALE.



Innovating. Educating. Changing Lives.™

New Standard of Care

Medical tubing doesn't have to slow you down. We've combined a soft and flexible silicone body with a breathable base material to give our patients comfort and mobility while holding their tubes securely in place.



Description	Order Number	Quantity
CINCH™	TSD0406-10	Box of 10
CINCH™	TSD0406-50	Box of 50
AMT Clamp™	4-3000	Box of 10

Device comes fully assembled with an easy peel away plastic backing

Innovating. Educating. Changing Lives.™

For a complete list of products visit us at www.AppliedMedical.net

800 869 7382

CS@AppliedMedical.net



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net

An ISO 13485 Registered Company CE