



Repositionable Clipping Device

A close-up photograph of a medical device. It consists of a long, black, flexible, corrugated tube that curves upwards. At the end of the tube is a blue cylindrical component, which is connected to a metal assembly. This assembly has four thin, metallic, curved arms or blades that are spread out in a star-like pattern. The background is a solid light blue.

Excellent Clamping & Rotation Performance
Shorter Implant Length

Patented product

Repositionable Clipping Device



CHARACTERISTICS

- *Patented Product*
- *Reposition Function*
- *Excellent Clamping Performance*
- *Excellent Rotation Performance*
- *Compatible MRI*
- *Greater Opening Width*
- *Shorter Implant Length*



Opening Width 16 mm	⇒	Implant Length 11.4 mm
Opening Width 12 mm	⇒	Implant Length 9.8 mm
Opening Width 9 mm	⇒	Implant Length 8.5 mm

SPECIFICATION

No.	Model	Working Length (mm)	Clip Opening (mm)	Jaw Angle (°)
1	WS-2727CD16135	2700	16	135
2	WS-2723CD16135	2300	16	135
3	WS-2719CD16135	1950	16	135
4	WS-2716CD16135	1650	16	135
5	WS-2727CD12135	2700	12	135
6	WS-2723CD12135	2300	12	135
7	WS-2719CD12135	1950	12	135
8	WS-2716CD12135	1650	12	135
9	WS-2727CD09135	2700	9	135
10	WS-2723CD09135	2300	9	135
11	WS-2719CD09135	1950	9	135
12	WS-2716CD09135	1650	9	135
13	WS-2727CD16090	2700	16	90
14	WS-2723CD16090	2300	16	90
15	WS-2719CD16090	1950	16	90
16	WS-2716CD16090	1650	16	90
17	WS-2727CD12090	2700	12	90
18	WS-2723CD12090	2300	12	90
19	WS-2719CD12090	1950	12	90
20	WS-2716CD12090	1650	12	90
21	WS-2727CD09090	2700	9	90
22	WS-2723CD09090	2300	9	90
23	WS-2719CD09090	1950	9	90
24	WS-2716CD09090	1650	9	90

* Max. Insertion Portion Diameter $\Phi 2.7\text{mm}$