

Disposable Injection Needles

Inject medicines into human tissue through endoscope.

Single use

CE 0197

Disposable Injection Needles



- Stainless steel needle point, Teflon needle inner tube and outer tube.
- Item No marked with '*' is our regular specification.



Item No.	Max. Insertion Portion Width & Working Length	Needle diameter	Needle length	Min. Insertion Channel Width	Sterile
WS-1816PN2304	Ø 1.9 × 1650 (mm)	23 G	4 mm	Ø 2.0 mm	●
WS-1816PN2504	Ø 1.9 × 1650 (mm)	25 G	4 mm	Ø 2.0 mm	●
WS-1816PN2305	Ø 1.9 × 1650 (mm)	23 G	5 mm	Ø 2.0 mm	●
WS-1816PN2505	Ø 1.9 × 1650 (mm)	25 G	5 mm	Ø 2.0 mm	●
WS-1816PN2306	Ø 1.9 × 1650 (mm)	23 G	6 mm	Ø 2.0 mm	●
WS-1816PN2308	Ø 1.9 × 1650 (mm)	23 G	8 mm	Ø 2.0 mm	●
WS-2416PN2104	Ø 2.5 × 1650 (mm)	21 G	4 mm	Ø 2.8 mm	●
WS-2416PN2304 *	Ø 2.5 × 1650 (mm)	23 G	4 mm	Ø 2.8 mm	●
WS-2416PN2504 *	Ø 2.5 × 1650 (mm)	25 G	4 mm	Ø 2.8 mm	●
WS-2416PN2105 *	Ø 2.5 × 1650 (mm)	21 G	5 mm	Ø 2.8 mm	●
WS-2416PN2505	Ø 2.5 × 1650 (mm)	25 G	5 mm	Ø 2.8 mm	●
WS-2416PN2306	Ø 2.5 × 1650 (mm)	23 G	6 mm	Ø 2.8 mm	●
WS-2416PN2506	Ø 2.5 × 1650 (mm)	25 G	6 mm	Ø 2.8 mm	●
WS-2416PN2108	Ø 2.5 × 1650 (mm)	21 G	8 mm	Ø 2.8 mm	●
WS-2423PN1904 *	Ø 2.5 × 2300 (mm)	19 G	4 mm	Ø 2.8 mm	●
WS-2423PN2304 *	Ø 2.5 × 2300 (mm)	23 G	4 mm	Ø 2.8 mm	●
WS-2423PN2504 *	Ø 2.5 × 2300 (mm)	25 G	4 mm	Ø 2.8 mm	●
WS-2423PN2105 *	Ø 2.5 × 2300 (mm)	21 G	5 mm	Ø 2.8 mm	●
WS-2423PN2305 *	Ø 2.5 × 2300 (mm)	23 G	5 mm	Ø 2.8 mm	●
WS-2423PN2505 *	Ø 2.5 × 2300 (mm)	25 G	5 mm	Ø 2.8 mm	●
WS-2423PN2306 *	Ø 2.5 × 2300 (mm)	23 G	6 mm	Ø 2.8 mm	●
WS-2423PN2506 *	Ø 2.5 × 2300 (mm)	25 G	6 mm	Ø 2.8 mm	●

