

# Disposable Electrosurgical Electrodes

Connect to high-frequency generator, cut and coagulate sphincter tissue.

Single use

## Disposable Needle Electrode • Triple lumen •



Cut tissue with high-frequency current\*\*




Item No.	Max. Insertion Portion Width & Working Length	Min. Insertion Channel Width	Needle length	Sterile
WS-2421DZ05	Ø 2.5×2100 (mm)	Ø 2.8 mm	5 mm	●
WS-2421DZ07	Ø 2.5×2100 (mm)	Ø 2.8 mm	7 mm	●
WS-2421DZ10	Ø 2.5×2100 (mm)	Ø 2.8 mm	10 mm	●

## Papillotome • One lumen •




Working Length: 2100 mm, Max. Insertion Portion Width: Ø2.5 mm, Min. Insertion Channel Width: Ø2.8 mm




Monofil Wire

Item No.	A	B	Sterile
WS-0020DJ1	0 mm	20 mm	●
WS-0025DJ1	0 mm	25 mm	●
WS-0030DJ1	0 mm	30 mm	●




Twisted Wire

Item No.	A	B	Sterile
WS-0020DJ1-B	0 mm	20 mm	●
WS-0025DJ1-B	0 mm	25 mm	●
WS-0030DJ1-B	0 mm	30 mm	●



Monofil Wire

Item No.	A	B	Sterile
WS-0020DJ1T	0 mm	20 mm	●
WS-0025DJ1T	0 mm	25 mm	●
WS-0030DJ1T	0 mm	30 mm	●



Twisted Wire

Item No.	A	B	Sterile
WS-0020DJ1T-B	0 mm	20 mm	●
WS-0025DJ1T-B	0 mm	25 mm	●
WS-0030DJ1T-B	0 mm	30 mm	●