



Disposable Grasping Forceps

Single use

3-Prong Type Grasping Forceps




Item No.	Max. Insertion Portion Width & Working Length	Min. Insertion Channel Width	Sterile
WS-2418GM3G	Ø 2.5 × 1800 (mm)	Ø 2.8 mm	●
WS-2423GM3G	Ø 2.5 × 2300 (mm)	Ø 2.8 mm	●




With Loop

Item No.	Max. Insertion Portion Width & Working Length	Min. Insertion Channel Width	Sterile
WS-2418GM3Q	Ø 2.5 × 1800 (mm)	Ø 2.8 mm	●
WS-2423GM3Q	Ø 2.5 × 2300 (mm)	Ø 2.8 mm	●

4-Prong Type Grasping Forceps



Item No.	Max. Insertion Portion Width & Working Length	Min. Insertion Channel Width	Sterile
WS-2418GM4G	Ø 2.5 × 1800 (mm)	Ø 2.8 mm	●
WS-2423GM4G	Ø 2.5 × 2300 (mm)	Ø 2.8 mm	●



With Loop

Item No.	Max. Insertion Portion Width & Working Length	Min. Insertion Channel Width	Sterile
WS-2418GM4Q	Ø 2.5 × 1800 (mm)	Ø 2.8 mm	●
WS-2423GM4Q	Ø 2.5 × 2300 (mm)	Ø 2.8 mm	●