

Disposable Endoscopic Distal Attachments



EN ISO: 13485/2012/AC:2012
CFDA: 辽械注准 20142220009

www.shangxianinc.com

Shangxian Minimal Invasive Inc.

Add:1st Floor, Block B2-2, China medicine innovation park,
Mulan road,Hi-tech development zone, Benxi 117004 China
Tel./Fax.: +86-24-6677-5326
E-mail: info@shangxianinc.com



Disposable Endoscopic Distal Attachments

Intended use

To be installed on the distal end of endoscope for maintaining good field of viewing.
Widely used in EMR, ESD, POEM surgeries and colonoscopy procedure.



Features

- Providing more working space for surgery.
- Maintaining the good field of viewing of endoscope procedure.
- when colonoscopy procedure, it helps to unfold the colon. More positive rate.
- Be aware of the distance between tissue and endoscope end.
- Tissue protection.
- Compatible to more endoscopes.



Specifications

Model	Distal Size	Working distance	Matching endoscope model
XT-DL-092-40	9.0~9.2mm	4.0mm	EG-2790K,EG-2770; GIF-Q260,GIF-H190
XT-DL-098-40	9.6~9.8mm	4.0mm	GIF-H260,GIF-H180,GIF-HQ190,PCF-PH190L/I
XT-DL-099-40	9.7~10.1mm	4.0mm	EC-2990Li,EC-2990Fi,EC-2990Mi,EG-2990K
XT-DL-105-40	10.3~10.7mm	4.0mm	PCF-Q260J,EC-3490TLi,EG-2990Zi,EG-2990K
XT-DL-108-40	10.6~11.0mm	4.0mm	GIF-H260Z,EG-590ZW,EG-450CT5,EG-2990i
XT-DL-115-40	11.3~11.7mm	4.0mm	PCF-H180AI/L,PCF-H190DI/L,EG-3470K, EC-3470MK/FK/LK,EG-3490K,EC-3490FK/Fi/LK,EG-450D5,
XT-DL-118-40	11.6~12.0mm	4.0mm	GIF-T2T40
XT-DL-122-40	12.0~12.4mm	4.0mm	CF-Q260A,CF-Q260D
XT-DL-128-40	12.6~13.0mm	4.0mm	EC-3890LZi/FZi/MZi,EC-380Kp,EC-250WM5, EC-590WM/ZW,EC-530WM/WI EC-450WM5/WI5
XT-DL-132-40	13.0~13.4mm	4.0mm	CF-Q180AL/I,CF-H190L/I,CF-HQ190L/I, EC-3890TFK/TLK,EC-3890Mi/Fi/Li,EG-3890TK
XT-DL-148-40	14.6~15.0mm	4.0mm	CF-H260AZL/I