

Disposable Hemostasis Clips & Delivery Systems

【Re-positionable】

- 15mm....=clip opening $\geq 15\text{m}$
- 360° sync-rotation exactly
- Proper clamping strength, no damage to tissue
- Sheath design protect endoscope working channel.




EN ISO: 13485/2012/AC:2012
CFDA: 辽械注准 20152220093

www.shangxianinc.com


Shangxian Minimal Invasive Inc.

Add:1st Floor, Block B2-2, China medicine innovation park,
Mulan road,Hi-tech development zone, Benxi 117004 China
Tel./Fax.: +86-24-6677-5326
E-mail: info@shangxianinc.com



Specifications

Models	Clip Opening	Working Length	Working Channel	Jaw Angle	Use For
SX-HR-1650-2	15mm	1650mm	≥2.8mm	90°±15°	胃镜
SX-HR-2300-2	15mm	2300mm	≥2.8mm	90°±15°	肠镜



Exact both sides rotation helps accurate clip deployment