


Disposable Hemostasis
Clips and Delivery
Systems(preload hemoclip)



EN ISO: 13485:2016
CFDA: 辽械注准20212020031

www.chinacleanday.cn

Shangxian Minimal Invasive Inc.
Add:1st Floor, Block B2-2, China medicine innovation park,
Mulan road,Hi-tech development zone, Benxi 117004 China
Tel: +86-24-6697 5946
Fax:+86-24-2378 2198

Disposable Hemostasis Clips and Delivery Systems

Intended use

Compatible with endoscope, intended for endoscopic clip placement within the gastrointestinal tract for the purpose of hemostasis for mucosal/submucosal defects, bleeding ulcers, arteries and polyps.

Features

- The clip is pre-load in disposable delivery system, which is convenient to use.
- Accurate sync-rotation offers exact clip deployment.
 - Reinforced design makes the clip performance more steady.
 - The biggest opening is 13/17mm after maximize the clip.
 - Multi times opening, 90/135 degree jaw.



Specifications

Model	Working length	Opening	Jaw	Working Channel
SX-HR-1650-2	1650mm	13mm	90°	≥2.8mm
SX-HR-2300-2	2300mm	13mm	90°	≥2.8mm
SX-HR-1650-1	1650mm	17mm	90°	≥2.8mm
SX-HR-2300-1	2300mm	17mm	90°	≥2.8mm