

Disposable Retrieval Net



Introduction

Used with endoscope, for retrieval of polyps and foreign bodies in the GI tract

Features

1. Large opening width
2. Slippery surface of the catheter glides nicely in the scope channel
3. Rotatable

Specifications

Model	Dia.(mm)	Working Length(mm)	Working Channel	Net diameter
MI-FG-TY-18*1600-15	1.8mm	1600mm	≥2.0mm	15mm
MI-FG-TY-18*1600-20	1.8mm	1600mm	≥2.0mm	20mm
MI-FG-TY-23*1600-20	2.3mm	1600mm	≥2.8mm	20mm
MI-FG-TY-23*1600-25	2.3mm	1600mm	≥2.8mm	25mm
MI-FG-TY-23*1600-30	2.3mm	1600mm	≥2.8mm	30mm
MI-FG-TY-23*1600-40	2.3mm	1600mm	≥2.8mm	40mm
MI-FG-TY-23*2300-20	2.3mm	2300mm	≥2.8mm	20mm
MI-FG-TY-23*2300-25	2.3mm	2300mm	≥2.8mm	25mm
MI-FG-TY-23*2300-30	2.3mm	2300mm	≥2.8mm	30mm
MI-FG-TY-23*2300-35	2.3mm	2300mm	≥2.8mm	35mm
MI-FG-TY-23*2300-40	2.3mm	2300mm	≥2.8mm	40mm