

Disposable Polypectomy Snares(hot & cold)



EN ISO: 13485:2016

www.chinacleanday.cn

Shangxian Minimal Invasive Inc.
Add:1st Floor, Block B2-2, China medicine innovation park,
Mulan road,Hi-tech development zone, Benxi 117004 China
Tel: +86-24-6697 5946
Fax:+86-24-2378 2198

Disposable Polypectomy Snares(hot & cold)



Intended use

For endoscopic clamp hole on the digestive tract polyps resection.

Features

- Good rotation performance and hand control feeling.
- Special alloyed material makes good shape retention.

Model	Dia.(mm)	Working channel(mm)	Working Length (mm)	Snare Diameter (mm)
SX-STY-18*2300-10(R)	1.8	≥2.0	2300	10
SX-STY-18*2300-15(R)	1.8	≥2.0	2300	15
SX-STY-18*2300-20(R)	1.8	≥2.0	2300	20
SX-STY-18*2300-25(R)	1.8	≥2.0	2300	25
SX-STY-18*2300-30(R)	1.8	≥2.0	2300	30
SX-STY-18*2300-35(R)	1.8	≥2.0	2300	35
SX-STY-18*2300-40(R)	1.8	≥2.0	2300	40
SX-STY-18*1600-10(R)	1.8	≥2.0	1600	10
SX-STY-18*1600-15(R)	1.8	≥2.0	1600	15
SX-STY-18*1600-20(R)	1.8	≥2.0	1600	20
SX-STY-18*1600-25(R)	1.8	≥2.0	1600	25
SX-STY-18*1600-30(R)	1.8	≥2.0	1600	30
SX-STY-18*1600-35(R)	1.8	≥2.0	1600	35
SX-STY-18*1600-40(R)	1.8	≥2.0	1600	40
SX-STY-23*2300-10(R)	2.3	≥2.8	2300	10
SX-STY-23*2300-15(R)	2.3	≥2.8	2300	15
SX-STY-23*2300-20(R)	2.3	≥2.8	2300	20
SX-STY-23*2300-25(R)	2.3	≥2.8	2300	25
SX-STY-23*2300-30(R)	2.3	≥2.8	2300	30
SX-STY-23*2300-35(R)	2.3	≥2.8	2300	35
SX-STY-23*2300-40(R)	2.3	≥2.8	2300	40
SX-STY-23*1600-10(R)	2.3	≥2.8	1600	10
SX-STY-23*1600-15(R)	2.3	≥2.8	1600	15
SX-STY-23*1600-20(R)	2.3	≥2.8	1600	20
SX-STY-23*1600-25(R)	2.3	≥2.8	1600	25
SX-STY-23*1600-30(R)	2.3	≥2.8	1600	30
SX-STY-23*1600-35(R)	2.3	≥2.8	1600	35
SX-STY-23*1600-40(R)	2.3	≥2.8	1600	40