

# Kiwi™ Omni-C Cup Vacuum Delivery System in C-Sections

The Kiwi™ Omni-C Cup is a complete vacuum delivery system designed specifically for use in confined space of the abdominal cavity.

The Kiwi Omni-C cup has approximately 30% lower profile than the Omnicup and a flexible stem, this allows easy insertion and manoeuvrability within the confined space of the abdominal cavity.

The two finger grooves for improved cup grip and also a baffle filter for vacuum dispersion facilitates proper placement on the flexion point.

The Kiwi Omni-C Cup is the ideal device for safe and effective removal of the baby through the uterine section with total control from the Kiwi PalmPump™.


With easy access to high floating head during repeat c-section, immediate assistance for failed vaginal deliveries when fetal head has descended low into the birth canal and, reduced pressure on the uterus.

Other clinical benefits from the Kiwi Omni-C Cup are reduced blood loss, smaller incisions, fewer hysterotomy extension and cervical lacerations. Without connections or tubing that could contaminate or interfere with the procedure.




### Ordering Information:

Part Number:	VAC-6000C
Quantity:	5 units per case
Delivery:	Usually ex-stock
Description:	Kiwi Omni-C Cup for C-Sections



Kiwi OmniCup on the left and the Kiwi Omni-C Cup on the right with profile approximately 30% lower.

Finger grooves and baffle filter



## KEBOMED

Formervangen 5, DK-2600 Glostrup  
(+45) 36398080, [info@kebomed.dk](mailto:info@kebomed.dk)  
[www.kebomed.dk](http://www.kebomed.dk)