


Disposable Endoscopic Hemoclip



Features

- Clip span width available in 9mm, 11mm, 13mm, 16mm, meet different clinical demands.
- Up to 5 times of repeated opening and closing prior to deployment, enable the repositioning of clip.
- 360° rotation, provide precise deployment.
- Unique jaws design provides strong clamping force.
- Short stem after deployment, create better visualization when multiple clips needed, and improve the comfort for the patient.

Product Code	Working Length (mm)	Diameter of Outer Tub (mm)	Max. Grasp Degree (°)	Max. Opening Size of Clip W (mm)	Endoscope Channel (mm)	Units/Box
CIT-ZXJ-(N)1650-9	1650	2.3	135°	9	≥ 2.8	10
CIT-ZXJ-(N)2350-9	2350	2.3	135°	9	≥ 2.8	10
CIT-ZXJ-(N)1650-11	1650	2.3	135°	11	≥ 2.8	10
CIT-ZXJ-(N)2350-11	2350	2.3	135°	11	≥ 2.8	10
CIT-ZXJ-(N)1650-13	1650	2.3	135°	13	≥ 2.8	10
CIT-ZXJ-(N)2350-13	2350	2.3	135°	13	≥ 2.8	10
CIT-ZXJ-(N)1650-16	1650	2.3	135°	16	≥ 2.8	10
CIT-ZXJ-(N)2350-16	2350	2.3	135°	16	≥ 2.8	10

Remark: code with (N) means without coating, code without (N) means with coating.