

GYNECARE TVT EXACT™ Continence System

Built upon many years of
GYNECARE TVT™ Success



GYNECARE TVT EXACT™

Tactile feedback for enhanced control

RIGID TROCAR

Rigid 3.0 mm diameter trocar shaft is designed to maintain control during passage while reducing penetration force.

SMOOTH, CLOSED-TIP TROCAR

Trocar features a smooth surface from handle to closed tip to reduce tissue drag during needle passage

ERGONOMIC HANDLE

Handle improves tactile feedback during tissue passage

Same dependable design

•The trocar curvature and tip radius of GYNECARE TVT™ Tension-Free Support for Incontinence and GYNECARE TVT EXACT™ Continence System are the same and designed so that the trocar maintains contact with posterior aspect of the pubic bone.¹⁴

•The Design of GYNECARE TVT™ with retropubic bottom to top approach has demonstrated significantly lower perforation and voiding dysfunction rate vs. top to bottom.¹²

•Laser-cut mesh.¹³



marketing@calderamedical.com

+1.818.879.6555



CALDERA
MEDICAL